



ATF-6 β Monoclonal Antibody

Catalog No	YP-mAb-01556
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ATF6B
Immunogen	The antiserum was produced against synthesized peptide derived from human ATF6B. AA range:401-450
Specificity	ATF-6 β Monoclonal Antibody detects endogenous levels of ATF-6 β protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Activating transcription factor 6 beta;ATF 6 beta;ATF6-beta;ATF6B; ATF6B_HUMAN;ATF6beta;cAMP response element binding protein related protein;cAMP response element-binding protein-related protein;cAMP responsive element binding protein like 1;cAMP-dependent transcription factor ATF-6 beta;cAMP-responsive element-binding protein-like 1;CREB L1;Creb related protein;CREB RP;Creb-rp;CREBL 1;CREBL1;Crebrp;Cyclic AMP dependent transcription factor ATF 6 beta;Cyclic AMP dependent transcription factor ATF6 beta;FLJ10066;G13;G13 protein;OTTHUMP00000029432; OTTHUMP00000029433;Processed cyclic AMP-dependent transcription factor ATF-6 beta;Protein G13;ATF-6 β ;Cyclic AMP-dependent transcription factor ATF-6 beta
Observed Band	90kD
Calculated Molecular Weight	90kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type II membrane protein . ; [Processed cyclic AMP-dependent transcription factor ATF-6 beta]: Nucleus . Under ER stress the cleaved N-terminal cytoplasmic domain translocates into the nucleus. .
Tissue Specificity	Ubiquitous.



Function

domain:The basic domain functions as a nuclear localization signal.,domain:The basic leucine-zipper domain is sufficient for association with the NF-Y trimer and binding to ERSE.,function:Transcriptional factor that acts in the unfolded protein response (UPR) pathway by activating UPR target genes induced during ER stress. Binds DNA on the 5'-CCAC[GA]-3' half of the ER stress response element (ERSE) (5'-CCAATN(9)CCAC[GA]-3') when NF-Y is bound to ERSE.,PTM:During unfolded protein response an approximative 60 kDa fragment containing the cytoplasmic transcription factor domain is released by proteolysis. The cleavage is probably performed sequentially by site-1 and site-2 proteases.,PTM:N-glycosylated.,similarity:Belongs to the bZIP family.,similarity:Belongs to the bZIP family. ATF subfamily.,similarity:Contains 1 bZIP domain.,subcellular location:Under ER stress the cleaved N-terminal cy

Background

The protein encoded by this gene is a transcription factor in the unfolded protein response (UPR) pathway during ER stress. Either as a homodimer or as a heterodimer with ATF6-alpha, the encoded protein binds to the ER stress response element, interacting with nuclear transcription factor Y to activate UPR target genes. The protein is normally found in the membrane of the endoplasmic reticulum; however, under ER stress, the N-terminal cytoplasmic domain is cleaved from the rest of the protein and translocates to the nucleus. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ATF-6β Monoclonal Antibody