



CAF-1 p60 Monoclonal Antibody

Catalog No	YP-mAb-01587
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHAF1B
Immunogen	The antiserum was produced against synthesized peptide derived from human CAF1B. AA range:71-120
Specificity	CAF-1 p60 Monoclonal Antibody detects endogenous levels of CAF-1 p60 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MPP 7;p60 CAF1;CAF1;CAF IP60;CAF I p60;CAF-1 p60;CHAF1B;Chromatin assembly factor 1 subunit B;CAF 1A;CAF1P60;MPHOSPH7;MPP7;CAF-1 subunit B;Chromatin assembly factor I p60 subunit;CAF-I 60 kDa subunit;CAF-I p60;M-phase phosphoprotein 7
Observed Band	61kD
Calculated Molecular Weight	61kD
Cell Pathway	Nucleus . Cytoplasm . DNA replication foci. Cytoplasmic in M phase.
Tissue Specificity	Epithelium,Lymphoblast,Uterus,
Function	developmental stage:Active complex is found in G1, S and G2 phases., function:Complex that is thought to mediate chromatin assembly in DNA replication and DNA repair. Assembles histone octamers onto replicating DNA in vitro. CAF-1 performs the first step of the nucleosome assembly process, bringing newly synthesized histones H3 and H4 to replicating DNA; histones H2A/H2B can bind to this chromatin precursor subsequent to DNA replication to complete the histone octamer. The CCR4-NOT complex functions as general transcription regulation complex.,PTM:Differentially phosphorylated during cell



cycle. During mitosis the p60 subunit of inactive CAF-1 is hyperphosphorylated and displaced into the cytosol. Progressively dephosphorylated from G1 to S and G2 phase. Phosphorylated p60 is recruited to chromatin undergoing DNA repair after UV irradiation in G1, S or G2 phases.,similarity:Belongs to the

Background

Chromatin assembly factor I (CAF-I) is required for the assembly of histone octamers onto newly-replicated DNA. CAF-I is composed of three protein subunits, p50, p60, and p150. The protein encoded by this gene corresponds to the p60 subunit and is required for chromatin assembly after replication. The encoded protein is differentially phosphorylated in a cell cycle-dependent manner. In addition, it is normally found in the nucleus except during mitosis, when it is released into the cytoplasm. This protein is a member of the WD-repeat HIR1 family and may also be involved in DNA repair. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CAF-1 p60 Monoclonal Antibody