



DNA pol γ Monoclonal Antibody

Catalog No	YP-mAb-01663
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POL1
Immunogen	The antiserum was produced against synthesized peptide derived from human POL1. AA range:641-690
Specificity	DNA pol γ Monoclonal Antibody detects endogenous levels of DNA pol γ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	DNA polymerase γ ;Eta2;POL1;polymerase (DNA directed) γ ;RAD30 homolog B;RAD30B;RAD30B;DNA pol γ
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus . Binding to ubiquitin mediates localization to replication forks after UV-induced DNA damage. .
Tissue Specificity	Ubiquitous. Highly expressed in testis.
Function	catalytic activity:Deoxynucleoside triphosphate + DNA(n) = diphosphate + DNA(n+1).,cofactor:Magnesium.,domain:The catalytic core consists of fingers, palm and thumb subdomains, but the fingers and thumb subdomains are much smaller than in high-fidelity polymerases; residues from five sequence motifs of the Y-family cluster around an active site cleft that can accommodate DNA and nucleotide substrates with relaxed geometric constraints, with consequently higher rates of misincorporation and low processivity.,function:Error-prone DNA polymerase specifically involved in DNA repair. Plays an important role in translesion synthesis, where the normal high-fidelity DNA polymerases cannot proceed and DNA synthesis stalls. Favors Hoogsteen base-pairing in the active



site. Inserts the correct base with high-fidelity opposite an adenosine template. Exhibits low fidelity and efficiency opposite a th

Background

catalytic activity: Deoxynucleoside triphosphate + DNA(n) = diphosphate + DNA(n+1)., cofactor: Magnesium., domain: The catalytic core consists of fingers, palm and thumb subdomains, but the fingers and thumb subdomains are much smaller than in high-fidelity polymerases; residues from five sequence motifs of the Y-family cluster around an active site cleft that can accommodate DNA and nucleotide substrates with relaxed geometric constraints, with consequently higher rates of misincorporation and low processivity., function: Error-prone DNA polymerase specifically involved in DNA repair. Plays an important role in translesion synthesis, where the normal high-fidelity DNA polymerases cannot proceed and DNA synthesis stalls. Favors Hoogsteen base-pairing in the active site. Inserts the correct base with high-fidelity opposite an adenosine template. Exhibits low fidelity and efficiency opposite a thymidine template, where it will preferentially insert guanosine. May play a role in hypermutation of immunoglobulin genes. Forms a Schiff base with 5'-deoxyribose phosphate at abasic sites, but may not have lyase activity., similarity: Belongs to the DNA polymerase type-Y family., similarity: Contains 1 umuC domain., subcellular location: Accumulates at replication forks after DNA damage., subunit: Binds REV1L (By similarity). Binds POLH., tissue specificity: Ubiquitous. Highly expressed in testis.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using DNA pol I Monoclonal Antibody