



FoxD4/D4L Monoclonal Antibody

Catalog No	YP-mAb-01722
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXD4
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXD4/L2/L3/L4/L5/L6. AA range:281-330
Specificity	FoxD4/D4L Monoclonal Antibody detects endogenous levels of FoxD4/D4L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FKHL9;forkhead box D4;Forkhead box protein D4;Forkhead related protein FKHL9;FOXD4;FOXD4A;FREAC 5;FREAC5;Myeloid factor alpha;FoxD4/D4L; Forkhead-related protein FKHL9;Forkhead-related transcription factor 5;FREAC-5;Myeloid factor-alpha;FOXD4L4;FOXD4B;Forkhead box protein D4-like 4; FOXD4-like 4;Forkhead box prote
Observed Band	46kD
Calculated Molecular Weight	46kD
Cell Pathway	Nucleus.
Tissue Specificity	Chondrosarcoma,Heart,PCR rescued clones,
Function	similarity:Contains 1 fork-head DNA-binding domain.,
Background	This gene encodes a member of the forkhead/winged helix-box (FOX) family of transcription factors. FOX transcription factors play critical roles in the regulation of multiple processes including metabolism, cell proliferation and gene expression during ontogenesis. Mutations in this gene are associated with a complex phenotype consisting of dilated cardiomyopathy, obsessive-compulsive



disorders, and suicidality. [provided by RefSeq, Mar 2012],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

