



GATA-5 Monoclonal Antibody

Catalog No	YP-mAb-01746
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GATA5
Immunogen	The antiserum was produced against synthesized peptide derived from human GATA-5. AA range:244-293
Specificity	GATA-5 Monoclonal Antibody detects endogenous levels of GATA-5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Transcription factor GATA-5;Transcription factor GATA 5;GATA-binding factor 5; GATA binding factor 5;GATA-5;GATA5
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Nucleus .
Tissue Specificity	Uterus,
Function	function:Binds to the functionally important CEF-1 nuclear protein binding site in the cardiac-specific slow/cardiac troponin C transcriptional enhancer. May play an important role in the transcriptional program(s) that underlies smooth muscle cell diversity.,similarity:Contains 2 GATA-type zinc fingers.,
Background	The protein encoded by this gene is a transcription factor that contains two GATA-type zinc fingers. The encoded protein is known to bind to hepatocyte nuclear factor-1alpha (HNF-1alpha), and this interaction is essential for cooperative activation of the intestinal lactase-phlorizin hydrolase promoter. In other organisms, similar proteins may be involved in the establishment of cardiac smooth muscle cell diversity. [provided by RefSeq, Jul 2008],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GATA-5 Monoclonal Antibody