



hnRNP CL1 Monoclonal Antibody

Catalog No	YP-mAb-01802
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HNRNPCL1
Immunogen	The antiserum was produced against synthesized peptide derived from human HNRCL. AA range:31-80
Specificity	hnRNP CL1 Monoclonal Antibody detects endogenous levels of hnRNP CL1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hnRNP CL1;HNRNPCL1;Heterogeneous nuclear ribonucleoprotein C-like 1; HNRPCL1;hnRNP C-like-1;hnRNP core protein C-like 1
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Nucleus . Component of ribonucleosomes. .
Tissue Specificity	Pooled,Uterus,
Function	function:May play a role in nucleosome assembly by neutralizing basic proteins such as A and B core hnRNPs.,similarity:Belongs to the RRM HNRPC family. RALY subfamily.,similarity:Contains 1 RRM (RNA recognition motif) domain., subcellular location:Component of ribonucleosomes.,
Background	function:May play a role in nucleosome assembly by neutralizing basic proteins such as A and B core hnRNPs.,similarity:Belongs to the RRM HNRPC family. RALY subfamily.,similarity:Contains 1 RRM (RNA recognition motif) domain., subcellular location:Component of ribonucleosomes.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

