



NLE1 Monoclonal Antibody

Catalog No	YP-mAb-01910
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NLE1
Immunogen	The antiserum was produced against synthesized peptide derived from human NLE1. AA range:301-350
Specificity	NLE1 Monoclonal Antibody detects endogenous levels of NLE1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Nle;NLE1;Notchless protein homolog 1;HUSSY-07
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Brain,Cervix carcinoma,Lung,Placenta,
Function	similarity:Contains 8 WD repeats.,
Background	similarity:Contains 8 WD repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

