



PRPF38A Monoclonal Antibody

Catalog No	YP-mAb-01962
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRPF38A
Immunogen	The antiserum was produced against synthesized peptide derived from human PRPF38A. AA range:104-153
Specificity	PRPF38A Monoclonal Antibody detects endogenous levels of PRPF38A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RP5 965L7.1;Pre-mRNA-splicing factor 38A;Pre mRNA splicing factor 38A; PRPF38A
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Cervix,Epithelium,Placenta,
Function	function:May be required for pre-mRNA splicing .,similarity:Belongs to the PRP38 family.,
Background	function:May be required for pre-mRNA splicing .,similarity:Belongs to the PRP38 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

