



RBM26 Monoclonal Antibody

Catalog No	YP-mAb-01969
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBM26
Immunogen	The antiserum was produced against synthesized peptide derived from human RBM26. AA range:711-760
Specificity	RBM26 Monoclonal Antibody detects endogenous levels of RBM26 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNA-binding motif protein 26;RBM 26;PRO1777;CTCL tumor antigen se70-2; C13orf10;RBM26;RNA-binding protein 26
Observed Band	113kD
Calculated Molecular Weight	113kD
Cell Pathway	
Tissue Specificity	Adipose tissue,Amygdala,Embryo,Epithelium,Fetal liver,Lymph,Placenta,Terato
Function	similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 2 RRM (RNA recognition motif) domains.,
Background	similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 2 RRM (RNA recognition motif) domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

