




Rrn3 Monoclonal Antibody

Catalog No	YP-mAb-01997
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RRN3
Immunogen	The antiserum was produced against synthesized peptide derived from human TIF-IA. AA range:602-651
Specificity	Rrn3 Monoclonal Antibody detects endogenous levels of Rrn3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNA polymerase I-specific transcription initiation factor RRN3;TIF IA;TIFIA;TIF-IA;RRN3;Transcription initiation factor IA
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Hippocampus,Kidney,Pooled,Spleen,
Function	function:Required for efficient transcription initiation by RNA polymerase I., similarity:Belongs to the RRN3 family.,
Background	function:Required for efficient transcription initiation by RNA polymerase I., similarity:Belongs to the RRN3 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Rrn3
Monoclonal Antibody