



Sox-6 Monoclonal Antibody

Catalog No	YP-mAb-02030
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SOX6
Immunogen	The antiserum was produced against synthesized peptide derived from human SOX6. AA range:10-59
Specificity	Sox-6 Monoclonal Antibody detects endogenous levels of Sox-6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSSOX6;SOXD;TOLCAS;Transcription factor SOX-6;Sox-6;SOX6
Observed Band	91kD
Calculated Molecular Weight	91kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Expressed in a wide variety of tissues, most abundantly in skeletal musclen.
Function	function:Transcriptional activator. Binds specifically to the DNA sequence 5'-AACAAAT-3'. Plays a key role in several developmental processes, including neurogenesis and skeleton formation.,PTM:Sumoylation inhibits the transcriptional activity.,similarity:Contains 1 HMG box DNA-binding domain., subunit:Interacts with DAZAP2.,tissue specificity:Expressed in a wide variety of tissues, most abundantly in skeletal muscle.,
Background	SRY-box 6(SOX6) Homo sapiens This gene encodes a member of the D subfamily of sex determining region y-related transcription factors that are characterized by a conserved DNA-binding domain termed the high mobility group box and by their ability to bind the minor groove of DNA. The encoded protein is a transcriptional activator that is required for normal development of



the central nervous system, chondrogenesis and maintenance of cardiac and skeletal muscle cells. The encoded protein interacts with other family members to cooperatively activate gene expression. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Mar 2009],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Sox-6 Monoclonal Antibody

