



SUHW3 Monoclonal Antibody


Catalog No	YP-mAb-02060
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF280C
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF280C. AA range:251-300
Specificity	SUHW3 Monoclonal Antibody detects endogenous levels of SUHW3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SUHW3;ZNF280C;Zinc finger protein 280C;ZNF633;Suppressor of hairy wing homolog 3;Zinc finger protein 633
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus .
Tissue Specificity	Colon,Testis,
Function	function:May function as a transcription factor.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 5 C2H2-type zinc fingers.,
Background	This gene encodes a member of the zinc finger domain-containing protein family. This family member contains multiple Cys2-His2(C2H2)-type zinc finger domains, the most common type of zinc finger domain that self-folds to form a beta-beta-alpha structure that binds a zinc ion. [provided by RefSeq, Aug 2011],
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SUHW3 Monoclonal Antibody