



TIF1α Monoclonal Antibody

Catalog No	YP-mAb-02116
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM24
Immunogen	The antiserum was produced against synthesized peptide derived from human TRIM24. AA range:1001-1050
Specificity	TIF1α Monoclonal Antibody detects endogenous levels of TIF1α protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TIF1-alpha;TIF1A;TIF1;RING-type E3 ubiquitin transferase TIF1-alpha; EC:2.3.2.27;TIF1α;TRIM24;Transcription intermediary factor 1-alpha;RNF82;E3 ubiquitin-protein ligase TRIM24;RING finger protein 82;Tripartite motif-containing protein 24
Observed Band	117kD
Calculated Molecular Weight	117kD
Cell Pathway	Nucleus . Cytoplasm . Mitochondrion . Colocalizes with sites of active transcription. Detected both in nucleus and cytoplasm in some breast cancer samples. Predominantly nuclear. Translocated from nucleus to mitochondria to mediate antiviral immunity (PubMed:32324863). .
Tissue Specificity	Brain,Epithelium,Mammary cancer,Testis,Thyroid,
Function	disease:A chromosomal aberration involving TIF1 is a cause of thyroid papillary carcinoma (PACT) [MIM:188550]. Translocation t(7;10)(q32;q11) with RET. The translocation generates the TIF1/RET (PTC6) oncogene.,function:Interacts selectively in vitro with the AF2-activating domain of the estrogen receptors. Association with DNA-bound estrogen receptors requires the presence of estradiol.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,



similarity:Contains 1 bromo domain.,similarity:Contains 1 PHD-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 2 B box-type zinc fingers.,subunit:Interacts with CBX1 and CBX3 (By similarity). Interacts with NR3C2.,

Background

The protein encoded by this gene mediates transcriptional control by interaction with the activation function 2 (AF2) region of several nuclear receptors, including the estrogen, retinoic acid, and vitamin D3 receptors. The protein localizes to nuclear bodies and is thought to associate with chromatin and heterochromatin-associated factors. The protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains - a RING, a B-box type 1 and a B-box type 2 - and a coiled-coil region. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using TIF1α Monoclonal Antibody