



Topo III β -1 Monoclonal Antibody

Catalog No	YP-mAb-02124
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TOP3B
Immunogen	The antiserum was produced against synthesized peptide derived from human TOP3B. AA range:51-100
Specificity	Topo III β -1 Monoclonal Antibody detects endogenous levels of Topo III β -1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Topo III β -1;TOP3B;DNA topoisomerase 3-beta-1;TOP3B1;DNA topoisomerase III beta-1
Observed Band	96kD
Calculated Molecular Weight	96kD
Cell Pathway	condensed chromosome,nucleus,cytosol,
Tissue Specificity	Isoform 1 is found in testis, heart and skeletal muscle. A 4 kb transcript which probably represents isoform 2 is found in thymus, kidney and pancreas.
Function	catalytic activity:ATP-independent breakage of single-stranded DNA, followed by passage and rejoining.,function:Possesses negatively supercoiled DNA relaxing activity.,similarity:Belongs to the prokaryotic type I/III topoisomerase family.,tissue specificity:Isoform 1 is found in testis, heart and skeletal muscle. A 4 kb transcript which probably represents isoform 2 is found in thymus, kidney and pancreas.,



Background

This gene encodes a DNA topoisomerase, an enzyme that controls and alters the topologic states of DNA during transcription. This enzyme catalyzes the transient breaking and rejoining of a single strand of DNA which allows the strands to pass through one another, thus relaxing the supercoils and altering the topology of DNA. The enzyme interacts with DNA helicase SGS1 and plays a role in DNA recombination, cellular aging and maintenance of genome stability. Low expression of this gene may be related to higher survival rates in breast cancer patients. This gene has a pseudogene on chromosome 22. Alternate splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Aug 2013],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Topo IIIβ-1 Monoclonal Antibody