



ZN-16 Monoclonal Antibody


Catalog No	YP-mAb-02176
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF292
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF292. AA range:1311-1360
Specificity	ZN-16 Monoclonal Antibody detects endogenous levels of ZN-16 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN-16;ZNF292;Zinc finger protein 292;KIAA0530
Observed Band	300kD
Calculated Molecular Weight	300kD
Cell Pathway	Nucleus .
Tissue Specificity	Bladder,Brain,Human endometrium,Placenta,Thalamus,
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZN-16 Monoclonal Antibody