



ZNF23 Monoclonal Antibody

Catalog No	YP-mAb-02179
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF23
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF23. AA range:151-200
Specificity	ZNF23 Monoclonal Antibody detects endogenous levels of ZNF23 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KOX16;ZNF359;ZNF612;Zfp612;ZNF23;Zinc finger protein 23;Zinc finger protein 359;Zinc finger protein 612;Zinc finger protein KOX16
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Nucleus .
Tissue Specificity	Atrial muscle,Brain,Hippocampus,Human small intestine,Lymphoid,Muscle, Substantia ni
Function	function:May be involved in transcriptional regulation. May have a role in embryonic development.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 17 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation. May have a role in embryonic development.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 17 C2H2-type zinc fingers.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

