



ZNF668 Monoclonal Antibody

Catalog No	YP-mAb-02198
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF668
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF668. AA range:271-320
Specificity	ZNF668 Monoclonal Antibody detects endogenous levels of ZNF668 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	zinc finger protein 668;ZNF668
Observed Band	68kD
Calculated Molecular Weight	68kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Mammary gland,Placenta,Skin,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 16 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 16 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZNF668 Monoclonal Antibody

