



ZSCAN16 Monoclonal Antibody


Catalog No	YP-mAb-02207
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZSCAN16
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF435. AA range:191-240
Specificity	ZSCAN16 Monoclonal Antibody detects endogenous levels of ZSCAN16 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZSCAN16;Zinc finger and SCAN domain-containing protein 16;ZNF392;ZNF435; Zinc finger protein 392;Zinc finger protein 435
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	Nucleus .
Tissue Specificity	Uterus,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 4 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 4 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZSCAN16 Monoclonal Antibody