



ZSCAN22 Monoclonal Antibody

Catalog No	YP-mAb-02208
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZSCAN22
Immunogen	The antiserum was produced against synthesized peptide derived from human ZSCAN22. AA range:321-370
Specificity	ZSCAN22 Monoclonal Antibody detects endogenous levels of ZSCAN22 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZSCAN22;Zinc finger and SCAN domain-containing protein 22;HKR2;ZNF50; Krueppel-related zinc finger protein 2;Protein HKR2;Zinc finger protein 50
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Nucleus .
Tissue Specificity	Fetal kidney,Placenta,Uterus,
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 8 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 8 C2H2-type zinc fingers.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

