



DNA Pol μ Monoclonal Antibody


Catalog No	YP-mAb-02227
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	POLM
Immunogen	Synthesized peptide derived from DNA Pol μ . at AA range: 210-290
Specificity	DNA Pol μ Monoclonal Antibody detects endogenous levels of DNA Pol μ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	DNA polymerase μ ;Pol Mu;POLM;polymerase (DNA directed); μ ;Tdt N;DNA Pol μ ;DNA-directed DNA/RNA polymerase μ ;polmu;Terminal transferase
Observed Band	54kD
Calculated Molecular Weight	54kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in a number of tissues. Abundant in thymus.
Function	catalytic activity:Deoxynucleoside triphosphate + DNA(n) = diphosphate + DNA(n+1).,cofactor:Magnesium.,function:Seems to act as an Ig mutase which is responsible for immunoglobulin (Ig) gene hypermutation.,similarity:Belongs to the DNA polymerase type-X family.,similarity:Contains 1 BRCT domain.,tissue specificity:Expressed in a number of tissues. Abundant in thymus.,
Background	catalytic activity:Deoxynucleoside triphosphate + DNA(n) = diphosphate + DNA(n+1).,cofactor:Magnesium.,function:Seems to act as an Ig mutase which is responsible for immunoglobulin (Ig) gene hypermutation.,similarity:Belongs to the DNA polymerase type-X family.,similarity:Contains 1 BRCT domain.,tissue specificity:Expressed in a number of tissues. Abundant in thymus.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using DNA Pol μ Monoclonal Antibody