



ABHD14A Monoclonal Antibody

Catalog No	YP-mAb-02453
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABHD14A
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD14A. AA range:201-250
Specificity	ABHD14A Monoclonal Antibody detects endogenous levels of ABHD14A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein ABHD14A;EC:3.-.-.;DORZ1;Alpha/beta hydrolase domain-containing protein 14A;Abhydrolase domain-containing protein 14A;ABHD14A
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Cytoplasm . Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Brain,Lymph,
Function	caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,
Background	caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

