



ABHD9 Monoclonal Antibody

Catalog No	YP-mAb-02460
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EPHX3
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD9. AA range:101-150
Specificity	ABHD9 Monoclonal Antibody detects endogenous levels of ABHD9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EH3;EC:3.3.2.1;Abhydrolase domain-containing protein 9;ABHD9;EPHX3; Epoxide hydrolase 3
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Microsome membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	EPHX3 (Epoxide Hydrolase 3) is a Protein Coding gene. GO annotations related to this gene include hydrolase activity. An important paralog of this gene is EPHX4.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

