



Atg16L2 Monoclonal Antibody

Catalog No	YP-mAb-02501
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ATG16L2
Immunogen	The antiserum was produced against synthesized peptide derived from human ATG16L2. AA range:241-290
Specificity	Atg16L2 Monoclonal Antibody detects endogenous levels of Atg16L2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	APG16 like 2;ATG16L2;Autophagy related protein 16 2;WDR80;Autophagy-related protein 16-2;APG16-like 2;WD repeat-containing protein 80
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Cytoplasm, cytosol . Localizes also to discrete punctae along the ciliary axoneme. .
Tissue Specificity	Astrocyte,Blood,Eye,Spleen,Stomach,
Function	function:May play a role in autophagy.,similarity:Belongs to the WD repeat ATG16 family.,similarity:Contains 7 WD repeats.,
Background	function:May play a role in autophagy.,similarity:Belongs to the WD repeat ATG16 family.,similarity:Contains 7 WD repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

