



CYP2D6 Monoclonal Antibody

Catalog No	YP-mAb-02586
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CYP2D6
Immunogen	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 2D6. AA range:251-300
Specificity	CYP2D6 Monoclonal Antibody detects endogenous levels of CYP2D6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CPD6;CYP2D;CYP2D@;CYP2D6;CYP2DL1;Cytochrome P450 2D6;P450 DB1;P450C2D;P450DB1;CYPIID6;Cytochrome P450-DB1;Debrisoquine 4-hydroxylase
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Endoplasmic reticulum membrane; Peripheral membrane protein. Microsome membrane ; Peripheral membrane protein.
Tissue Specificity	Brain,Liver,PCR rescued clones,
Function	catalytic activity:RH + reduced flavoprotein + O(2) = ROH + oxidized flavoprotein + H(2)O.,cofactor:Heme group.,function:Responsible for the metabolism of many drugs and environmental chemicals that it oxidizes. It is involved in the metabolism of drugs such as antiarrhythmics, adrenoceptor antagonists, and tricyclic antidepressants.,induction:By pregnancy.,online information:CYP2D6 alleles,online information:CYP2D6 entry,polymorphism:Allele CYP2D6*7 was also known as CYP2D6E, allele CYP2D6*9 as CYP2D6C, allele CYP2D6*10 as CYP2D6J, allele CYP2D6*17 as CYP2D6Z.,polymorphism:Genetic variations in CYP2D6 are the cause of poor drug metabolism CYP2D6-related [MIM:608902].



The CYP2D6 gene is highly polymorphic. CYP2D6 activity ranges widely within a population comprising ultrarapid (UM), extensive (EM), intermediate (IM) and poor (PM) metabolizer phenotypes. UM and PM are those most at risk f

Background

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and is known to metabolize as many as 25% of commonly prescribed drugs. Its substrates include antidepressants, antipsychotics, analgesics and antitussives, beta adrenergic blocking agents, antiarrhythmics and antiemetics. The gene is highly polymorphic in the human population; certain alleles result in the poor metabolizer phenotype, characterized by a decreased ability to metabolize the enzyme's substrates. Some individuals with the poor metabolizer phenotype have no functional protein since they carry 2 null alleles whereas in other individuals the gene is absent. This gene can vary in


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CYP2D6 Monoclonal Antibody