



DGAT2L6 Monoclonal Antibody

Catalog No	YP-mAb-02615
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DGAT2L6
Immunogen	The antiserum was produced against synthesized peptide derived from human DGAT2L6. AA range:271-320
Specificity	DGAT2L6 Monoclonal Antibody detects endogenous levels of DGAT2L6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DGAT2L6;Diacylglycerol O-acyltransferase 2-like protein 6;DC3;Diacylglycerol O-acyltransferase candidate 3;hDC3
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Endoplasmic reticulum membrane; Multi-pass membrane protein .
Tissue Specificity	Expressed in all tissues tested except pancreas.
Function	function:Probable acyltransferase uses fatty acyl-CoA as substrate (By similarity) . Has no wax synthase activity to produce wax esters.,similarity:Belongs to the diacylglycerol acyltransferase family.,tissue specificity:Expressed in all tissues tested except pancreas.,
Background	This gene is a member of the diacylglycerol acyltransferase 2 family. The encoded protein is a putative acyltransferase and is most likely involved in the synthesis of di- or triacylglycerol, however its substrate specificity is currently unknown. [provided by RefSeq, Feb 2015],



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

