



ELOVL5 Monoclonal Antibody

Catalog No	YP-mAb-02630
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ELOVL5
Immunogen	The antiserum was produced against synthesized peptide derived from human ELOVL5. AA range:250-299
Specificity	ELOVL5 Monoclonal Antibody detects endogenous levels of ELOVL5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ELOVL2;ELOVL5;Fatty acid elongase 5;HELO1;Elongation of very long chain fatty acids protein 5;3-keto acyl-CoA synthase ELOVL5;ELOVL fatty acid elongase 5;ELOVL FA elongase 5;Fatty acid elongase 1
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Cell projection, dendrite . In Purkinje cells, the protein localizes to the soma and proximal portion of the dendritic tree. .
Tissue Specificity	Ubiquitous. Highly expressed in the adrenal gland and testis. Weakly expressed in prostate, lung and brain. Expressed in the cerebellum.
Function	function:May be involved in the elongation of long-chain polyunsaturated fatty acids (PUFAs).,similarity:Belongs to the ELO family.,tissue specificity:Highly expressed in the adrenal gland and testis. Weekly expressed in prostate, lung and brain.,
Background	This gene belongs to the ELO family. It is highly expressed in the adrenal gland and testis, and encodes a multi-pass membrane protein that is localized in the endoplasmic reticulum. This protein is involved in the elongation of long-chain



polyunsaturated fatty acids. Mutations in this gene have been associated with spinocerebellar ataxia-38 (SCA38). Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2014],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

