



FAM35A/B Monoclonal Antibody

Catalog No	YP-mAb-02635
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM35A/FAM35B
Immunogen	The antiserum was produced against synthesized peptide derived from human FAM35A/B. AA range:394-443
Specificity	FAM35A/B Monoclonal Antibody detects endogenous levels of FAM35A/B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM35A/B;FAM35A/FAM35B;Protein FAM35A/B;FAM35A;Protein FAM35A;FAM35B;Protein FAM35B
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	Nucleus . Chromosome . Localizes to nuclear foci in response to DNA damage. .
Tissue Specificity	Brain,Eye,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the FAM35 family.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the FAM35 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

