



Mfn2 Monoclonal Antibody

Catalog No	YP-mAb-02675
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MFN2
Immunogen	The antiserum was produced against synthesized peptide derived from human Mfn2. AA range:354-403
Specificity	Mfn2 Monoclonal Antibody detects endogenous levels of Mfn2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CPRP1;EC:3.6.5.-;HSG;KIAA0214;MARF;Mfn2;Mitofusin-2;Transmembrane GTPase MFN2
Observed Band	86kD
Calculated Molecular Weight	86kD
Cell Pathway	Mitochondrion outer membrane ; Multi-pass membrane protein . Colocalizes with BAX during apoptosis. .
Tissue Specificity	Ubiquitous; expressed at low level. Highly expressed in heart and kidney.
Function	catalytic activity:GTP + H(2)O = GDP + phosphate.,disease:Defects in MFN2 are the cause of Charcot-Marie-Tooth disease type 2A2 (CMT2A2) [MIM:609260] . CMT2A2 is a form of Charcot-Marie-Tooth disease, the most common inherited disorder of the peripheral nervous system. Charcot-Marie-Tooth disease is classified in two main groups on the basis of electrophysiologic properties and histopathology: primary peripheral demyelinating neuropathy or CMT1, and primary peripheral axonal neuropathy or CMT2. Neuropathies of the CMT2 group are characterized by signs of axonal regeneration in the absence of obvious myelin alterations, normal or slightly reduced nerve conduction velocities, and progressive distal muscle weakness and atrophy.,



disease: Defects in MFN2 are the cause of Charcot-Marie-Tooth disease type 6 (CMT6) [MIM:601152]; also referred to as autosomal dominant hereditary motor and sensory n

Background

This gene encodes a mitochondrial membrane protein that participates in mitochondrial fusion and contributes to the maintenance and operation of the mitochondrial network. This protein is involved in the regulation of vascular smooth muscle cell proliferation, and it may play a role in the pathophysiology of obesity. Mutations in this gene cause Charcot-Marie-Tooth disease type 2A2, and hereditary motor and sensory neuropathy VI, which are both disorders of the peripheral nervous system. Defects in this gene have also been associated with early-onset stroke. Two transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Mfn2 Monoclonal Antibody