



NQO1 Monoclonal Antibody

Catalog No	YP-mAb-02723
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NQO1
Immunogen	The antiserum was produced against synthesized peptide derived from human NQO1. AA range:203-252
Specificity	NQO1 Monoclonal Antibody detects endogenous levels of NQO1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Azoreductase;DHQU;DIA 4;DIA4;DT diaphorase;NQO1;NAD(P)H dehydrogenase [quinone] 1;NMOR1;DT-diaphorase;DTD;Menadione reductase; NAD(P)H:quinone oxidoreductase 1;Phylloquinone reductase;Quinone reductase 1;QR1
Observed Band	31kD
Calculated Molecular Weight	31kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Colon,Liver,Pooled,
Function	catalytic activity:NAD(P)H + a quinone = NAD(P)(+) + a hydroquinone., cofactor:FAD.,enzyme regulation:Inhibited by dicoumarol.,function:The enzyme apparently serves as a quinone reductase in connection with conjugation reactions of hydroquinons involved in detoxification pathways as well as in biosynthetic processes such as the vitamin K-dependent gamma-carboxylation of glutamate residues in prothrombin synthesis.,induction:By dioxin.,mass spectrometry: PubMed:11735396,miscellaneous:Quinone reductase accepts electrons from both NADH and NADPH with equal efficiency.,polymorphism:The Ser-187 polymorphism may be linked to susceptibility to forms of cancers.,



similarity:Belongs to the NAD(P)H dehydrogenase (quinone) family.,
subunit:Homodimer.,

Background

This gene is a member of the NAD(P)H dehydrogenase (quinone) family and encodes a cytoplasmic 2-electron reductase. This FAD-binding protein forms homodimers and reduces quinones to hydroquinones. This protein's enzymatic activity prevents the one electron reduction of quinones that results in the production of radical species. Mutations in this gene have been associated with tardive dyskinesia (TD), an increased risk of hematotoxicity after exposure to benzene, and susceptibility to various forms of cancer. Altered expression of this protein has been seen in many tumors and is also associated with Alzheimer's disease (AD). Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using NQO1 Monoclonal Antibody