



p44S10 Monoclonal Antibody

Catalog No	YP-mAb-02730
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSMD6
Immunogen	The antiserum was produced against synthesized peptide derived from human PSMD6. AA range:201-250
Specificity	p44S10 Monoclonal Antibody detects endogenous levels of p44S10 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0107;p42A;p44S10;PFAAP4;PSMD6;Rpn7;S10;SGA 113M;26S proteasome non-ATPase regulatory subunit 6;26S proteasome regulatory subunit RPN7;26S proteasome regulatory subunit S10;Breast cancer-associated protein SGA-113M;Phosphonoformate immuno-associated protein 4
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	proteasome complex,nucleoplasm,cytosol,proteasome regulatory particle, proteasome accessory complex,extracellular exosome,
Tissue Specificity	Bone marrow,Lymph,Mammary cancer,Ovary,Placenta,Testis,
Function	function:Acts as a regulatory subunit of the 26S proteasome which is involved in the ATP-dependent degradation of ubiquitinated proteins.,similarity:Belongs to the proteasome subunit S10 family.,similarity:Contains 1 PCI domain., subunit:Component of the PA700 complex.,
Background	This gene encodes a member of the protease subunit S10 family. The encoded protein is a subunit of the 26S proteasome which colocalizes with DNA damage foci and is involved in the ATP-dependent degradation of ubiquinated proteins.



Alternative splicing results in multiple transcript variants [provided by RefSeq, Nov 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using p44S10 Monoclonal Antibody