




UBPH Monoclonal Antibody

Catalog No	YP-mAb-02818
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UBFD1
Immunogen	The antiserum was produced against synthesized peptide derived from human UBFD1. AA range:135-184
Specificity	UBPH Monoclonal Antibody detects endogenous levels of UBPH protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ubiquitin domain-containing protein UBFD1;Ubiquitin Family Domain Containing 1;Ubiquitin-binding protein homolog;UBPH;UBFD1
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	cell-cell adherens junction,
Tissue Specificity	Skeletal muscle,
Function	
Background	UBFD1 (Ubiquitin Family Domain Containing 1) is a Protein Coding gene. GO annotations related to this gene include poly (A) RNA binding.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using UBPH Monoclonal Antibody