



USP13 Monoclonal Antibody

Catalog No	YP-mAb-02827
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP13
Immunogen	The antiserum was produced against synthesized peptide derived from human USP13. AA range:811-860
Specificity	USP13 Monoclonal Antibody detects endogenous levels of USP13 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.4.19.12;Isopeptidase T-3;ISOT3;ISOT-3;Ubiquitin carboxyl-terminal hydrolase 13;USP13;Deubiquitinating enzyme 13;Ubiquitin thioesterase 13; Ubiquitin-specific-processing protease 13
Observed Band	97kD
Calculated Molecular Weight	97kD
Cell Pathway	
Tissue Specificity	Highly expressed in ovary and testes.
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBP-type zinc finger.,similarity:Contains 2 UBA domains.,tissue specificity:Highly expressed in ovary and testes.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBP-type zinc finger.,similarity:Contains 2 UBA domains.,tissue specificity:Highly expressed in ovary and testes.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using USP13 Monoclonal Antibody