



Actin α 3 Monoclonal Antibody

Catalog No	YP-mAb-03078
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACTG2
Immunogen	The antiserum was produced against synthesized peptide derived from human Actin-gamma2. AA range:1-50
Specificity	Actin α 3 Monoclonal Antibody detects endogenous levels of Actin α 3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Actin binding protein;Actinin alpha 2;ACTN 2;ACTN2;ACTN2_HUMAN;Alpha actinin 2;Alpha actinin skeletal muscle;Alpha actinin skeletal muscle isoform 2; Alpha-actinin skeletal muscle isoform 2;Alpha-actinin-2;CMD1AA;F actin cross linking protein;F-actin cross-linking protein;Actinin alpha 3;actinin;alpha;ACTN 3; ACTN3;ACTN3_HUMAN;Alpha actinin skeletal muscle isoform 3;Alpha-Actinin 3; Alpha-actinin skeletal muscle isoform 3;Alpha-actinin-3;MGC1172;MGC1175; OTTHUMP21949;OTTHUMP21954;Actin α 3;ACTG2;Actin gamma-enteric smooth muscle;ACTA3;ACTL3;ACTSG;Actin;gamma-enteric smooth muscle; Alpha-actin-3;Gamma-2-actin;Smooth muscle gamma-actin
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Cytoplasm, cytoskeleton.
Tissue Specificity	In the intestine, abundantly expressed in smooth muscle cells of muscularis mucosa and muscularis propria. Also detected in intestinal vascular smooth muscle cells.



Function

function: Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells., miscellaneous: In vertebrates 3 main groups of actin isoforms, alpha, beta and gamma have been identified. The alpha actins are found in muscle tissues and are a major constituent of the contractile apparatus. The beta and gamma actins coexist in most cell types as components of the cytoskeleton and as mediators of internal cell motility., similarity: Belongs to the actin family., subunit: Polymerization of globular actin (G-actin) leads to a structural filament (F-actin) in the form of a two-stranded helix. Each actin can bind to 4 others.,

Background

Actins are highly conserved proteins that are involved in various types of cell motility and in the maintenance of the cytoskeleton. Three types of actins, alpha, beta and gamma, have been identified in vertebrates. Alpha actins are found in muscle tissues and are a major constituent of the contractile apparatus. The beta and gamma actins co-exist in most cell types as components of the cytoskeleton and as mediators of internal cell motility. This gene encodes actin gamma 2; a smooth muscle actin found in enteric tissues. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Based on similarity to peptide cleavage of related actins, the mature protein of this gene is formed by removal of two N-terminal peptides.[provided by RefSeq, Dec 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Actin α 3 Monoclonal Antibody