




Moesin/Ezrin/Radixin Monoclonal Antibody

Catalog No	YP-mAb-03152
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MSN
Immunogen	The antiserum was produced against synthesized peptide derived from human Moesin/Ezrin/Radixin. AA range:524-573
Specificity	Moesin/Ezrin/Radixin Monoclonal Antibody detects endogenous levels of Moesin/Ezrin/Radixin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Villin 2 ezrin;CVIL;CVL;Cytovillin 2;Cytovillin;DKFZp762H157;Epididymis secretory protein Li 15;EZR;EZRI_HUMAN;Ezrin;FLJ26216;HEL S 15;MGC1584; p81;VIL 2;VIL2;Villin 2 (ezrin);Villin 2;Villin-2;Villin2;Epididymis luminal protein 7; HEL7;Membrane organizing extension spike protein;Membrane-organizing extension spike protein;MOES_HUMAN;Moesin;Moesin/anaplastic lymphoma kinase fusion protein;Msn;MSN/ALK fusion;CB567;CG12537;DFNB24;ESP1;Hh-induced MATH and BTB domain-containing protein;HIB;Moesin-B;Protein roadkill;RADI_HUMAN;Radixin;RDX;Moesin/Ezrin/Radixin
Observed Band	67kD
Calculated Molecular Weight	67kD
Cell Pathway	Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cytoplasm, cytoskeleton . Apical cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell projection, microvillus membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell projection, microvillus . Phosphorylated form is enriched in microvilli-like structures at apical membrane. Increased cell membrane localization of both phosphorylated and non-phosphorylated forms seen after thrombin treatment (By similarity). Localizes at the uropods of T lymphoblasts. .



Tissue Specificity	In all tissues and cultured cells studied.
Function	function:Probably involved in connections of major cytoskeletal structures to the plasma membrane.,PTM:Phosphorylation on Thr-558 is crucial for the formation of microvilli-like structures.,similarity:Contains 1 FERM domain.,subcellular location:Phosphorylated form is enriched in microvilli-like structures at apical membrane.,subunit:In resting T-cells, part of a PAG1-SLC9A3R1-MSN complex which is disrupted upon TCR activation (By similarity). Binds SLC9A3R1.,tissue specificity:In all tissues and cultured cells studied.,
Background	Moesin (for membrane-organizing extension spike protein) is a member of the ERM family which includes ezrin and radixin. ERM proteins appear to function as cross-linkers between plasma membranes and actin-based cytoskeletons. Moesin is localized to filopodia and other membranous protrusions that are important for cell-cell recognition and signaling and for cell movement. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Moesin/Ezrin/Radixin Monoclonal Antibody