



TCP-1 η Monoclonal Antibody

Catalog No	YP-mAb-03195
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCT7
Immunogen	The antiserum was produced against synthesized peptide derived from human CCT7. AA range:221-270
Specificity	TCP-1 η Monoclonal Antibody detects endogenous levels of TCP-1 η protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	TCP-1 η ;CCT7;T-complex protein 1 subunit eta;CCTH;NIP7-1;TCP-1-eta;CCT-eta;HIV-1 Nef-interacting protein
Observed Band	59kD
Calculated Molecular Weight	59kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Adipose tissue,Brain,Cajal-Retzius cell,Cerebellum,Colon carcinoma,
Function	function:Molecular chaperone; assist the folding of proteins upon ATP hydrolysis. Known to play a role, in vitro, in the folding of actin and tubulin., similarity:Belongs to the chaperonin (HSP60) family.,similarity:Belongs to the TCP-1 chaperonin family.,subunit:Heterooligomeric complex of about 850 to 900 kDa that forms two stacked rings, 12 to 16 nm in diameter. Interacts with PACRG.,
Background	This gene encodes a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex



and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 5 and 6. [provided by RefSeq, Oct 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using TCP-1 η Monoclonal Antibody