



mPR α Monoclonal Antibody


Catalog No	YP-mAb-03316
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PAQR7
Immunogen	The antiserum was produced against synthesized peptide derived from human MPRA. AA range:297-346
Specificity	mPR α Monoclonal Antibody detects endogenous levels of mPR α protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	mPR α ;PAQR7;Membrane progestin receptor alpha;MRPA;mPR alpha;Progestin and adipoQ receptor family member 7;Progestin and adipoQ receptor family member VII
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in a wide range of tissues including ovary, testis, placenta, uterus and bladder.
Function	function:Steroid membrane receptor. Binds progesterone in vitro. May be involved in oocyte maturation.,similarity:Belongs to the ADIPOR family.,tissue specificity:Expressed in ovary, testis, placenta, uterus and bladder.,
Background	function:Steroid membrane receptor. Binds progesterone in vitro. May be involved in oocyte maturation.,similarity:Belongs to the ADIPOR family.,tissue specificity:Expressed in ovary, testis, placenta, uterus and bladder.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using mPR α Monoclonal Antibody