



RXR γ Monoclonal Antibody

Catalog No	YP-mAb-03335
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RXRG
Immunogen	The antiserum was produced against synthesized peptide derived from human Retinoid X Receptor gamma. AA range:171-220
Specificity	RXR γ Monoclonal Antibody detects endogenous levels of RXR γ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ28;FLJ318;FLJ162;FLJ16733;MGC1272;NR2B1;Nuclear receptor subfamily 2 group B member 1;OTTHUMP2251;Retinoic acid receptor RXR alpha;Retinoic acid receptor RXR-alpha;Retinoid X nuclear receptor alpha;Retinoid X receptor alpha;RXR alpha1;Rxra;RXRA_HUMAN;RXRalpha1;RXR γ ;RXRG;Retinoic acid receptor RXR-gamma;NR2B3;Nuclear receptor subfamily 2 group B member 3; Retinoid X receptor gamma
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Expressed in aortic endothelial cells (at protein level).
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., domain:Composed of three domains: a modulating N-terminal domain, a DNA-binding domain and a C-terminal steroid-binding domain.,function:Nuclear hormone receptor. Involved in the retinoic acid response pathway. Binds 9-cis retinoic acid (9C-RA).,similarity:Belongs to the nuclear hormone receptor family. NR2 subfamily.,similarity:Contains 1 nuclear receptor DNA-binding domain.,



Background

retinoid X receptor gamma(RXRγ) Homo sapiens This gene encodes a member of the retinoid X receptor (RXR) family of nuclear receptors which are involved in mediating the antiproliferative effects of retinoic acid (RA). This receptor forms dimers with the retinoic acid, thyroid hormone, and vitamin D receptors, increasing both DNA binding and transcriptional function on their respective response elements. This gene is expressed at significantly lower levels in non-small cell lung cancer cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jun 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using RXRγ Monoclonal Antibody