




Arrdc4 Monoclonal Antibody

Catalog No	YP-mAb-03721
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARRDC4
Immunogen	The antiserum was produced against synthesized peptide derived from human ARRDC4. AA range:127-176
Specificity	Arrdc4 Monoclonal Antibody detects endogenous levels of Arrdc4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARRDC4_HUMAN;ARRDC 4;Arrdc4;Arrestin domain containing protein 4;Arrestin domain-containing protein 4;FLJ3645
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Early endosome . Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cytoplasmic vesicle . Also found in extracellular vesicles different from exosomes. .
Tissue Specificity	Placenta,Testis,Thalamus,
Function	similarity:Belongs to the arrestin family.,
Background	similarity:Belongs to the arrestin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Arrdc4 Monoclonal Antibody