



FGF-BP2 Monoclonal Antibody


Catalog No	YP-mAb-03871
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FGFBP2
Immunogen	The antiserum was produced against synthesized peptide derived from human FGFP2. AA range:170-219
Specificity	FGF-BP2 Monoclonal Antibody detects endogenous levels of FGF-BP2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HBp17-related protein;Fibroblast growth factor-binding protein 2;FGFBP-2;FGF-BP2;FGF-binding protein 2;FGFBP2;KSP37;37 kDa killer-specific secretory protein;HBp17-RP
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Secreted, extracellular space .
Tissue Specificity	Expressed in serum, peripheral leukocytes and cytotoxic T-lymphocytes, but not in granulocytes and monocytes (at protein level).
Function	similarity:Belongs to the fibroblast growth factor-binding protein family.,tissue specificity:Expressed in serum, peripheral leukocytes and cytotoxic T lymphocytes, but not in granulocytes and monocytes (at protein level).,
Background	This gene encodes a member of the fibroblast growth factor binding protein family. The encoded protein is a serum protein that is selectively secreted by cytotoxic lymphocytes and may be involved in cytotoxic lymphocyte-mediated immunity. An increase in the amount of gene product may be associated with atopic asthma and mild extrinsic asthma.[provided by RefSeq Staff, Oct 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FGF-BP2 Monoclonal Antibody