



HSP β 2 Monoclonal Antibody

Catalog No	YP-mAb-03918
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HSPB2
Immunogen	The antiserum was produced against synthesized peptide derived from human HSPB2. AA range:121-170
Specificity	HSP β 2 Monoclonal Antibody detects endogenous levels of HSP β 2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	HSP β 2;HSPB2;Heat shock protein beta-2;DMPK-binding protein;MKBP
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Cytoplasm . Nucleus . Localizes to nuclear foci.
Tissue Specificity	Expressed preferentially in skeletal muscle and heart but not in the lens.
Function	similarity:Belongs to the small heat shock protein (HSP20) family., subunit:Associates specifically with myotonic dystrophy protein kinase (DMPK)., tissue specificity:Expressed preferentially in skeletal muscle and heart but not in the lens.,
Background	The protein encoded by this gene belongs to the superfamily of small heat-shock proteins containing a conservative alpha-crystallin domain at the C-terminal part of the molecule. The protein is expressed preferentially in the heart and skeletal muscle. This protein regulates Myotonic Dystrophy Protein Kinase, which plays an important role in maintenance of muscle structure and function. [provided by RefSeq, Dec 2012],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using HSP β 2 Monoclonal Antibody