



MRP7 Monoclonal Antibody

Catalog No	YP-mAb-03972
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABCC10
Immunogen	The antiserum was produced against synthesized peptide derived from human MRP7. AA range:835-884
Specificity	MRP7 Monoclonal Antibody detects endogenous levels of MRP7 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABCC7;ATP-binding cassette sub-family C member 7;cAMP-dependent chloride channel;CFTR/MRP;Channel conductance-controlling ATPase;MRP7;ABCC10; Multidrug resistance-associated protein 7;SIMRP7;ATP-binding cassette sub-family C member 10
Observed Band	162kD
Calculated Molecular Weight	162kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Basolateral cell membrane ; Multi-pass membrane protein .
Tissue Specificity	[Isoform 1]: Specifically expressed in spleen. ; [Isoform 2]: Widely expressed.
Function	function:ATP-dependent transporter probably involved in cellular detoxification through lipophilic anion extrusion.,similarity:Belongs to the ABC transporter family.,similarity:Belongs to the ABC transporter family. Conjugate transporter (TC 3.A.1.208) subfamily.,similarity:Contains 2 ABC transmembrane type-1 domains.,similarity:Contains 2 ABC transporter domains.,tissue specificity:Isoform 1 is specifically expressed in spleen. Isoform 2 is more widely expressed.,



Background

The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, and White). This ABC full-transporter is a member of the MRP subfamily which is involved in multi-drug resistance. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Nov 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP7 Monoclonal Antibody