



MRP-S18C Monoclonal Antibody

Catalog No	YP-mAb-04008
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPS18C
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPS18C. AA range:71-120
Specificity	MRP-S18C Monoclonal Antibody detects endogenous levels of MRP-S18C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	28S ribosomal protein S18-1;mitochondrial;BS18m;CGI 134;CGI134;CGI-134;MRP-S18C;MRPS18C;28S ribosomal protein S18c mitochondrial;28S ribosomal protein S18c;MRP-S18-c;Mrps18-c;S18mt-c;MRP-S18-1; mitochondrial
Observed Band	16kD
Calculated Molecular Weight	16kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Brain,
Function	similarity:Belongs to the ribosomal protein S18P family..subunit:Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian



mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that belongs to the ribosomal protein S18P family. The encoded protein is one of three that has significant sequence similarity to bacterial S18 proteins. The primary sequences of the three human mitochondria

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

