



Rabr Monoclonal Antibody

Catalog No	YP-mAb-04108
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AGFG2
Immunogen	The antiserum was produced against synthesized peptide derived from human AGFG2. AA range:91-140
Specificity	Rabr Monoclonal Antibody detects endogenous levels of Rabr protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RABR;RAB R;HRBL;HIV-1 Rev-binding protein-like protein;ArfGAP with FG repeats 2;AGFG2;Arf-GAP domain and FG repeat-containing protein 2;Rev/Rex activation domain-binding protein related;RAB-R
Observed Band	49kD
Calculated Molecular Weight	49kD
Cell Pathway	membrane,
Tissue Specificity	Brain,Kidney,Placenta,
Function	domain:Contains FG repeats and 4 N-P-F repeats.,similarity:Contains 1 Arf-GAP domain.,subunit:Interacts with EPS15R.,
Background	ArfGAP with FG repeats 2(AGFG2) Homo sapiens This gene is a member of the HIV-1 Rev binding protein (HRB) family and encodes a protein with one Arf-GAP zinc finger domain, several phe-gly (FG) motifs, and four asn-pro-phe (NPF) motifs. This protein interacts with Eps15 homology (EH) domains and plays a role in the Rev export pathway, which mediates the nucleocytoplasmic transfer of proteins and RNAs. Alternatively spliced transcript variants have been described, but their biological validity has not been determined. [provided by



RefSeq, Feb 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

