



RB3 Monoclonal Antibody

Catalog No	YP-mAb-04111
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STMN4
Immunogen	The antiserum was produced against synthesized peptide derived from human STMN4. AA range:140-189
Specificity	RB3 Monoclonal Antibody detects endogenous levels of RB3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RB3;Stathmin 4;stathmin like 4;Stathmin like protein B3;STMN4;Stathmin-4;Stathmin-like protein B3
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Golgi apparatus . Cell projection, growth cone . Cell projection, axon .
Tissue Specificity	Amygdala,Brain,
Function	similarity:Belongs to the stathmin family.,
Background	similarity:Belongs to the stathmin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

