



RBQ-3 Monoclonal Antibody

Catalog No	YP-mAb-04112
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBBP5
Immunogen	Synthesized peptide derived from RBQ-3 . at AA range: 190-270
Specificity	RBQ-3 Monoclonal Antibody detects endogenous levels of RBQ-3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RBBP 5;RBBP5;RBQ3;SWD1;RBQ-3;Retinoblastoma-binding protein 5;RBBP-5; Retinoblastoma-binding protein RBQ-3
Observed Band	59kD
Calculated Molecular Weight	59kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	function:Binds preferentially to underphosphorylated retinoblastoma protein., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 6 WD repeats.,subunit:Component of the SET1 complex, at least composed of the catalytic subunit (SETD1A or SETD1B), WDR5, WDR82, RBBP5, ASH2/ASH2L and CXXC1/CFP1. Component of MLL-containing complexes (named MLL, ASCOM, MLL2/MLL3 or MLL3/MLL4 complex): at least composed ASH2L, RBBP5, DPY30, WDR5, one or several histone methyltransferases (MLL, MLL2, MLL3 and/or MLL4), and the facultative components MEN1, HCFC1, HCFC2, NCOA6, KDM6A, PAXIP1/PTIP and C16orf53/PA1.,tissue specificity:Ubiquitously expressed.,



Background

This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. The encoded protein binds directly to retinoblastoma protein, which regulates cell proliferation. It interacts preferentially with the underphosphorylated retinoblastoma protein via the E1A-binding pocket B. Three alternatively spliced transcript variants that encode different protein isoforms have been described for this gene. [provided by RefSeq, Jul 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using RBQ-3 Monoclonal Antibody