



Ribosomal Protein L23 Monoclonal Antibody

Catalog No	YP-mAb-04134
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL23
Immunogen	The antiserum was produced against synthesized peptide derived from human RPL23. AA range:51-100
Specificity	Ribosomal Protein L23 Monoclonal Antibody detects endogenous levels of Ribosomal Protein L23 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Large ribosomal subunit protein uL14;ribosomal protein L23;rpL17;RPL23;60S ribosomal protein L23;60S ribosomal protein L17
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	nucleolus,cytoplasm,cytosol,ribosome,focal adhesion,membrane,cytosolic large ribosomal subunit,extracellular matrix,extracellular exosome,
Tissue Specificity	Brain,Kidney,Uterus,
Function	similarity:Belongs to the ribosomal protein L14P family.,
Background	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L14P family of ribosomal proteins. It is located in the cytoplasm. This gene has been referred to as rpL17 because the encoded protein shares amino acid identity with ribosomal protein L17 from <i>Saccharomyces cerevisiae</i> ;



however, its official symbol is RPL23. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Ribosomal Protein L23 Monoclonal Antibody