



Ribosomal Protein S15 Monoclonal Antibody

Catalog No	YP-mAb-04157
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS15
Immunogen	The antiserum was produced against synthesized peptide derived from human RPS15. AA range:41-90
Specificity	Ribosomal Protein S15 Monoclonal Antibody detects endogenous levels of Ribosomal Protein S15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small ribosomal subunit protein uS19;RIG protein;RIG;ribosomal protein S15; RPS15;40S ribosomal protein S15
Observed Band	17kD
Calculated Molecular Weight	17kD
Cell Pathway	nucleus,nucleoplasm,cytoplasm,cytosol,ribosome,focal adhesion,small ribosomal subunit,membrane,cytosolic small ribosomal subunit,
Tissue Specificity	Pancreas,Placenta,
Function	similarity:Belongs to the ribosomal protein S19P family.,
Background	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S19P family of ribosomal proteins. It is located in the cytoplasm. This gene has been found to be activated in various tumors, such as insulinomas, esophageal cancers, and colon cancers. As is typical for genes



encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Ribosomal Protein S15 Monoclonal Antibody