




Serpin A11 Monoclonal Antibody

Catalog No	YP-mAb-04203
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SERPINA11
Immunogen	The antiserum was produced against synthesized peptide derived from human SERPINA11. AA range:191-240
Specificity	Serpin A11 Monoclonal Antibody detects endogenous levels of Serpin A11 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Serpin A11;antiproteinase-like 2;SERPINA11
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	Secreted .
Tissue Specificity	Fetal liver,
Function	similarity:Belongs to the serpin family.,
Background	similarity:Belongs to the serpin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Serpin A11 Monoclonal Antibody