



SH3BGR Monoclonal Antibody

Catalog No	YP-mAb-04206
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SH3BGR
Immunogen	The antiserum was produced against synthesized peptide derived from human SH3BGR. AA range:81-130
Specificity	SH3BGR Monoclonal Antibody detects endogenous levels of SH3BGR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	21 GARP;21 glutamic acid rich protein;SH3BGR;SH3BGR protein;SH3 domain-binding glutamic acid-rich protein;21-glutamic acid-rich protein;21-GARP
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	cytosol,
Tissue Specificity	Expressed in heart and skeletal muscle.
Function	caution:It is uncertain whether Met-1 or Met-64 is the initiator.,similarity:Belongs to the SH3BGR family.,tissue specificity:Expressed in heart and skeletal muscle.,
Background	caution:It is uncertain whether Met-1 or Met-64 is the initiator.,similarity:Belongs to the SH3BGR family.,tissue specificity:Expressed in heart and skeletal muscle.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using SH3BGR Monoclonal Antibody

