



SVOP Monoclonal Antibody

Catalog No	YP-mAb-04230
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SVOP
Immunogen	The antiserum was produced against synthesized peptide derived from human SVOP. AA range:271-320
Specificity	SVOP Monoclonal Antibody detects endogenous levels of SVOP protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLC22B4;SVOP;Synaptic vesicle 2-related protein;SV2-related protein
Observed Band	61kD
Calculated Molecular Weight	61kD
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane ; Multi-pass membrane protein .
Tissue Specificity	Amygdala,Brain,
Function	similarity:Belongs to the major facilitator superfamily.,
Background	similarity:Belongs to the major facilitator superfamily.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

